

# Supercritical Fluid Application Notes

**SCF  
511**

## Extraction of PAHs from Soil Using Supercritical Fluids

[Back to List](#)

### Introduction

Supercritical CO<sub>2</sub> alone easily dissolves polyaromatic hydrocarbons (PAHs) from spiked matrices, but not from weathered, aged samples. The addition of a small amount of co-solvent is necessary to displace the PAHs from the surface of a weathered sample for complete extraction.

The addition of a fixed quantity of methanol to standard reference material prior to supercritical CO<sub>2</sub> extraction shows increased recoveries of PAHs when compared to US EPA liquid/solid extraction techniques.

### Equipment

- ✓ Applied Separations' *Spe-ed* SFE Supercritical Extraction System



### Materials

- ✓ Standard Reference Material (PAHs in Soil)
- ✓ SFE grade CO<sub>2</sub>
- ✓ C18 Collection cartridge (1g/6mL #12007)
- ✓ Pesticide grade methanol

### Method

Place 0.5g of soil/sediment into a 1mL extraction vessel and add 100µL of methanol. Extract sample at specified conditions.

### Extraction Conditions

Extraction vessel:	1mL
Sample:	0.5 g
Pressure:	7000 psi
Temperature:	80°C
CO <sub>2</sub> Flow Rate:	2L/min (gas)
Collection:	SPE cartridge 1g/6mL C18
Static:	10 minutes
Dynamic:	30 minutes
SPE Rinse:	5mL methanol spiked with internal standard

### Results

SFE Recoveries of PAHs from Soil Certified by Soxhlet

	Certified (mg/kg)	SFE (mg/kg)
Naphthalene	32	31
Acenaphthene	19	17
Phenanthrene	1618	1573
Anthracene	422	493
Benzo(bk)flouranthene	152	180

**Applied  
Separations**

Distributor:



ERC GmbH

☎ +49 89 66055696  
Fax +49 89 60824826  
partner@erc-hplc.de  
www.erc-hplc.de



# Supercritical Fluid Application Notes

## Analysis

GC-MS

## Conclusion

A fixed quantity of co-solvent addition to a soil sample prior to extraction gave excellent results with 100 $\mu$ L of methanol added. Only small quantities of methanol are required to displace the PAH analytes from the soil matrix and do not desorb the analyte from the solid phase extraction cartridge used to trap the analyte.

---

Applied  
Separations

Distributor:



ERC GmbH

+49 89 66055696

Fax +49 89 60824826

partner@erc-hplc.de

www.erc-hplc.de