

# Antec Zubehör und Ersatzteile ROXY™ und DECADE II

| VT-03 Cell  |                |  |
|-------------|----------------|--|
| Bestell Nr. | Hersteller Nr. | Zellen mit GlassyCarbon Arbeitselektrode |
| 03005       | 110.4105       | VT-03 cell 3 mm GC sb                    |
| 03011       | 110.4205       | VT-03 cell 3 mm GC ISAAC                 |
| 03017       | 110.4305       | VT-03 cell 3 mm GC HyREF                 |
| 03006       | 110.4115       | VT-03 cell 2 mm GC sb                    |
| 03012       | 110.4215       | VT-03 cell 2 mm GC ISAAC                 |
| 03018       | 110.4315       | VT-03 cell 2 mm GC HyREF                 |
| 03023       | 111.4105       | VT-03 cell 0.7 mm GC sb                  |
| 03024       | 111.4205       | VT-03 cell 0.7 mm GC ISAAC               |
| 03025       | 111.4305       | VT-03 cell 0.7 mm GC HyREF               |
| Bestell Nr. | Hersteller Nr. | Zellen mit Metallarbeitselektrode        |
| 03007       | 110.4120       | VT-03 cell 3 mm Pt sb                    |
| 03013       | 110.4220       | VT-03 cell 3 mm Pt ISAAC                 |
| 03019       | 110.4320       | VT-03 cell 3 mm Pt HyREF                 |
| 03008       | 110.4125       | VT-03 cell 3 mm Au sb                    |
| 03014       | 110.4225       | VT-03 cell 3 mm Au ISAAC                 |
| 03020       | 110.4325       | VT-03 cell 3 mm Au HyREF                 |
| 03009       | 110.4130       | VT-03 cell 2 mm Ag sb                    |
| 03015       | 110.4230       | VT-03 cell 2 mm Ag ISAAC                 |
| 03021       | 110.4330       | VT-03 cell 2 mm Ag HyREF                 |
| 03010       | 110.4135       | VT-03 cell 3 mm Cu sb                    |
| 03016       | 110.4235       | VT-03 cell 3 mm Cu ISAAC                 |
| 03022       | 110.4335       | VT-03 cell 3 mm Cu HyREF                 |
| Bestell Nr. | Hersteller Nr. | Zubehör                                  |
| 03079       | 110.2011       | Spacer for VT-03 cell, 120 µm            |
| 03080       | 110.2010       | Spacer for VT-03 cell, 50 µm             |
| 03081       | 110.2009       | Spacer for VT-03 cell, 25 µm             |
| 03026       | 110.1300       | sb REF complete                          |
| 03091       | 250.1025       | Polishing disc for WE                    |
| 03095       | 250.1026       | Polishing disc for REF                   |
| 03097       | 250.1040       | Flattening plate 12 µm                   |
| 03098       | 250.1042       | Flattening plate 30 µm                   |
| 03099       | 250.1045       | flattening/polishing kit for metal WE    |
| 03092       | 250.1030       | 10 mL diamond slurry 1 µm                |
| 03100       | 250.2005       | 50ml KCl solution sat'd, with AgCl       |
| 03101       | 250.2010       | ISAAC solution 10 ml                     |
| 03082       | 250.1575       | Fingertight fitting VT-03, M5            |

# Antec **Zubehör und Ersatzteile ROXY™ und DECADE II**

| Flexcell    |                |                                 |
|-------------|----------------|---------------------------------|
| Bestell Nr. | Hersteller Nr. |                                 |
| 03032       | 102.4305M      | Flexcell MD HyREF               |
| 03029       | 102.4305       | Flexcell GC HyREF               |
| 03030       | 102.4320       | Flexcell Pt HyREF               |
| 03027       | 102.4325       | Flexcell Au HyREF               |
| 03031       | 102.4330       | Flexcell Ag HyREF               |
|             |                |                                 |
| 03033       | 102.4105       | Flexcell GC sb                  |
| 03034       | 102.4120       | Flexcell Pt sb                  |
| 03035       | 102.4125       | Flexcell Au sb                  |
| 03036       | 102.4130       | Flexcell Ag sb                  |
|             |                |                                 |
| 03037       | 102.4205       | Flexcell GC ISAAC               |
| 03038       | 102.4220       | Flexcell Pt ISAAC               |
| 03039       | 102.4225       | Flexcell Au ISAAC               |
| 03040       | 102.4230       | Flexcell Ag ISAAC               |
| Bestell Nr. | Hersteller Nr. | <b>Arbeitselektroden</b>        |
| 03051       | 102.5007       | WE disc GC                      |
| 03052       | 102.5022       | WE disc Pt                      |
| 03053       | 102.5027       | WE disc Au                      |
| 03054       | 102.5032       | WE disc Ag                      |
| 03055       | 102.5037       | WE disc Cu                      |
| 03056       | 210.5050       | WE disc MD                      |
| Bestell Nr. | Hersteller Nr. | <b>Zubehör</b>                  |
| 03078       | 102.2217       | Spacer for F/R cell, 50 µm      |
| 03102       | 102.0913       | HyREF for F/R                   |
| 03091       | 250.1025       | Polishing disc for WE           |
| 03092       | 250.1030       | 10 mL diamond slurry 1 µm       |
| 03083       | 250.1571       | Fingertight fitting PCTFE 10-32 |
| 03026       | 110.1300       | sb REF complete                 |
| 03103       | 102.0914       | ISAAC for F/R                   |

# Antec **Zubehör und Ersatzteile ROXY™ und DECADE II**

| Reactor Cell                        |                |                                     |
|-------------------------------------|----------------|-------------------------------------|
| Bestell Nr.                         | Hersteller Nr. |                                     |
| 03104                               | 210.0040       | Reactor cell kit (incl MD,GC,Au,Pt) |
|                                     |                |                                     |
| Bestell Nr.                         | Hersteller Nr. | <b>Arbeitselektroden</b>            |
| 03051                               | 210.5007       | WE disc GC                          |
| 03052                               | 210.5022       | WE disc Pt                          |
| 03053                               | 210.5027       | WE disc Au                          |
| 03054                               | 210.5032       | WE disc Ag                          |
| 03056                               | 210.5050       | WE disc MD                          |
| Bestell Nr.                         | Hersteller Nr. | <b>Zubehör</b>                      |
| 03078                               | 210.2217       | Spacer 50 µm for Reactor            |
| 03102                               | 210.0913       | HyREF for ReactorCell               |
| 03091                               | 250.1025       | Polishing disc for WE               |
| 03092                               | 250.1030       | 10 mL diamond slurry 1 µm           |
| 03083                               | 250.1571       | Fingertight fitting PCTFE 10-32     |
| µPrep Cell                          |                |                                     |
| Bestell Nr.                         | Hersteller Nr. |                                     |
| 03063                               | 204.4300       | µ-PrepCell incl. Kit                |
| Bestell Nr.                         | Hersteller Nr. | <b>Arbeitselektroden</b>            |
| 03105                               | 204.5007       | µ-PrepCell WE GC                    |
| 03106                               | 204.5050       | µ-PrepCell WE MD                    |
| Bestell Nr.                         | Hersteller Nr. | <b>Zubehör</b>                      |
| 03107                               | 204.0500       | µ-PrepCell O-ring inlet block       |
| 03108                               | 204.0502       | µ-PrepCell O-ring REF electrode     |
| 03109                               | 204.2217       | µ-PrepCell spacer 50 µm             |
| 03110                               | 204.2218       | µ-PrepCell spacer 100 µm            |
| 03091                               | 250.1025       | Polishing disc for WE               |
| 03092                               | 250.1030       | 10 mL diamond slurry 1 µm           |
| 03097                               | 250.1040       | Flattening plate 12 µm              |
| 03098                               | 250.1042       | Flattening plate 30 µm              |
| 03111                               | 250.1570       | Fingertight fitting PEEK 10-32      |
| <b>Weiteres Zubehör auf Anfrage</b> |                |                                     |