

SPE

Solid Phase
Extraction

Application Note

World Leaders in Sample Preparation

SPE105

ISOLATION OF AMINE ANTIOXIDANT FROM MOTOR OIL

Spe-ed™ Cartridge	Cat. No. 2107- Silica Gel, 1000mg/6mL. Process with <i>Spe-ed Mate</i>
Sample Preparation	Mix 1gm. oil and 1mL hexane into a vial.
Sample Addition	Add prepared sample to cartridge. Wash sample vial with 1mL of hexane and add to cartridge. Aspirate sample through cartridge.
Cartridge Wash	6mL hexane Discard elute.
Analyte Elution	4mL methylene chloride



**Applied
Separations**

Distributor:



ERC GmbH
☎ +49 89 66055696
Fax +49 89 60824826
partner@erc-hplc.de
www.erc-hplc.de