

# Treiber für Sedex ELSD's



		SEDEX 85	SEDEX 90	SEDEX 100	SEDEX FP	SEDEX LC
Anbieter	Software					
Agilent	ChemStation Edition™	✓	✓	✓	✓	✓
	OpenLab ChemStation Edition™	✓	✓	✓	✓	✓
	EZChrom Elite™	✓	✓	✓	✓	✓
	OpenLab EZChrom Edition™	✓	✓	✓	✓	✓
DataApex	Clarity™	✓	✓	✓	✓	✓
Knauer	PurityChrom™5				✓	✓
	ChromGate™	✓	✓			
Thermo Fisher (Dionex)	Chromeleon™6.8	✓	✓			
	Chromeleon™7.2	✓	✓	✓	✓	✓
	ChromQuest™	✓	✓			
	Xcalibur™	✓	✓			
SCPA	PrepCon™5				✓	✓
	Chromstar™7				✓	✓